

SPECIAL STANDARD BORE HOSE ASSEMBLY

Hose Assy Oxygen, 1.5m, 9/16" UNF (DISS) Nut to (DISS) Nut

PHOTOGRAPH REQUIRED

PRODUCT CODE

AHD02**SPECIFICATION**

- British Standard end fittings (Where Applicable)
- British Standard Colour Coded Hose (Length advised in product code)
- Anti Static and Phthalate Free

OVERVIEW

Oxylitre Holdings Limited offer the most comprehensive range of Pipeline Hose Assemblies. Our Pipeline Hose Assemblies Range cover all Medical Gases from oxygen through to special gas mixtures. Made from flexible anti-static phalate free colour coordinated hose. Oxylitre also offer a full design service, to create the hose assembly you require, for more information on this service please contact sales. Our listed hose assemblies comprise of an extensive range of gas specific connections. A full and comprehensive look at our range can be found in our product index.

Each Pipeline Hose Assembly is produced to BS and European Standards and all are available in none listed lengths. Oxylitre's Hose Assembly Range comes with a minimum 5 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom
+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

