

AMBULANCE PIPELINE SYSTEMS

Line Pressure Audible & Visual Warning System

PRODUCT CODE

APW100



✓ 4 bar

✓ 90dB

✓ 12Volt

SPECIFICATION

- Fire retardant & IP65 rated ergonomic enclosure
- Nominal 4 bar input pressure
- Audio Alarm at 90dB
- Pipeline Status LED's
- 3/8" BSP Male Connection for Pipeline Nut and Tail
- 12Volt DC
- Low Pressure Output Signal

OVERVIEW

The Oxylitre' Line Pressure Warning System, have been designed to work alongside any ambulance pipeline system, specifically designed to give an audible and visual warning of low pressure in the pipeline system. The system includes an on and off switch also pipeline status LED's the system comes complete with a direct feed from the low pressure signal which can be connected to an Ambulance Control Module.

Oxylitre's Line Pressure Warning System is dependable, easily maintained with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

