

STANDARD MEDICAL PIPELINE FLOWMETER

Standard Single Flowmeter Medical Air 4 Bar 0-15lpm Complete with Humidifier Bottle

PRODUCT CODE

F1611H



✓ 4 Bar

✓ 0-15Lpm

SPECIFICATION

- Medical Air 4 Bar (MA-4)
- Brass Chrome Plated Body
- NOT MRI Conditional
- 9/16 UNF Fitting
- British Standard Probe
- 4 bar nominal pressure
- BS Colour Fine adjustment Control
- Humidifier Bottle (Autoclavable to 121°C) complete with Tubing Nipple
- 0-15Lpm (MA-4)

WEIGHTS AND DIMENSIONS

Weight	0.38 kg
Height	305 mm
Width	76 mm
Depth	149 mm

ACCESSORIES

F1611H

SP46 Tubing Nipple female 9/16 UNF (White for O2)

H110 Humidifier (Bubble Type) Autoclavable 250ml Bottle

IS1174 Oxygen Connecting Tube (1.8m)

SKF6-02 5 Year Service Pack

OVERVIEW

Designed around the needs of Healthcare and Emergency Service Providers, Oxylitre's range of Standard Medical Pipeline Flowmeters are the Flowmeter of choice for delivering flexible solutions. Our Standard Medical Pipeline Flowmeter Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Medical Pipeline Flowmeters offer high performance with durability engineered in to the robust design.

Oxylitre's range of Standard Medical Pipeline Flowmeters are dependable, extremely easy to maintain with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

