

STANDARD MEDICAL PIPELINE FLOWMETER

Standard Single Flowmeter Rail Mounted Oxygen 0- 15lpm with Humidifier Bottle

PRODUCT CODE

F1630-H



✓ 4 bar

✓ 0-15Lpm

SPECIFICATION

- OXYGEN
- Brass Chrome Plated Body
- NOT MRI Conditional
- 9/16 UNF Fitting
- 4 bar nominal pressure
- BS Colour Fine adjustment Control
- Humidifier Bottle (Autoclavable to 121°C) complete with Tubing Nipple
- 0-15Lpm (O₂)
- Complete with Universal Rail Clamp
- 1.5 Metre BS Coloured Hose and Probe

WEIGHTS AND DIMENSIONS

Weight	0.77 kg
Height	300 mm
Width	44 mm
Depth	148 mm

ACCESSORIES

F1630-H

- SP46** Tubing Nipple female 9/16 UNF (White for Oxygen)
- H110** Humidifier (Bubble Type) Autoclavable 250ml Bottle
- SKF6-O2** 5 Year Service Pack

OVERVIEW

Designed around the needs of Healthcare and Emergency Service Providers, Oxylitre's range of Standard Medical Pipeline Flowmeters are the Flowmeter of choice for delivering flexible solutions. Our Standard Medical Pipeline Flowmeter Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Medical Pipeline Flowmeters offer high performance with durability engineered in to the robust design.

Oxylitre's range of Standard Medical Pipeline Flowmeters are dependable, extremely easy to maintain with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

