

STANDARD MEDICAL PIPELINE FLOWMETER

Standard Twin Flowmeters Medical Air 4 Bar 0-15Lpm

PRODUCT CODE

F1632A



✓ 4 Bar

✓ 0-15Lpm

SPECIFICATION

- Medical Air 4 Bar (MA-4)
- Brass Chrome Plated Body
- NOT MRI Conditional
- 9/16 UNF Fitting
- 4 bar nominal pressure
- BS Colour Fine adjustment Control
- Air Safety Shield & Tubing Nipple
- 0-15Lpm (MA-4)
- "Y" Configuration Twin Mounting Block to British Standard Probe
- Complete with Universal Rail Clamp
- 1.5 Metre BS Coloured Hose and Probe

WEIGHTS AND DIMENSIONS

Weight	1.15 kg
Height	204 mm
Width	208 mm
Depth	155 mm

ACCESSORIES

F1632A

SP46 Tubing Nipple female 9/16 UNF (White for Oxygen)

H110 Humidifier (Bubble Type) Autoclavable 250ml Bottle

SKF6-O2 5 Year Service Pack

OVERVIEW

Designed around the needs of Healthcare and Emergency Service Providers, Oxylitre's range of Standard Medical Pipeline Flowmeters are the Flowmeter of choice for delivering flexible solutions. Our Standard Medical Pipeline Flowmeter Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Medical Pipeline Flowmeters offer high performance with durability engineered in to the robust design.

Oxylitre's range of Standard Medical Pipeline Flowmeters are dependable, extremely easy to maintain with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 · sales@oxylitre.co.uk · service@oxylitre.co.uk · oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

