

MEDICAL SELECT-A-FLOW PIPELINE FLOWMETER

Select-A-Flow Flowmeter (Oxygen) 0-15Lpm and Self Sealing Valve

PRODUCT CODE

F16705



SPECIFICATION

- Select-A-Flow Flowmeter
- 4 bar nominal pressure
- BS Colour coded "Click Style" control
- 11 Individual Flow Settings
- Tubing Nipple
- Standard Back Bar Mounting to British Standard Probe
- Complete with Single BS Coloured Self Sealing Valve

ACCESSORIES

F16705

IS1174 Connecting Tube (1.8m)

SKF8-15LPM 5 Year Service Pack

OVERVIEW

Oxylitre's range of Medical Select-A-Flow (SAF) Pipeline Flowmeters are the Flowmeter of choice for delivering flexible solutions for Emergency Services and Healthcare Providers. Our Medical (SAF) Pipeline Flowmeter Range are simple to use, reliable and functional with 11 flow settings as standard. The Medical (SAF) Pipeline Flowmeters offer high performance with durability engineered in to the robust and compact design.

Oxylitre's range of Medical Select-A-Flow Pipeline Flowmeters are dependable, extremely easy to maintain with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom
+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

