

COMPACT MEDICAL PIPELINE FLOWMETER

Compact Single Flowmeter Medical Air 4 Bar 0-15Lpm 3/8" Nut & Liner Connection

PRODUCT CODE

F2600A



✓ 4 Bar

✓ 0-15Lpm

SPECIFICATION

- Medical Air 4 Bar (MA-4)
- Moulded BS Colour Coded Body
- MRI Conditional
- 9/16 UNF Fitting
- 3/8" Nut & Liner Connection
- 4 bar nominal pressure
- BS Colour Fine adjustment Control
- Air Safety Shield & Tubing Nipple
- 0-15Lpm (MA-4)

ACCESSORIES

F2600A

SP46 Tubing Nipple female 9/16 UNF (White for O2)

H110 Humidifier (Bubble Type) Autoclavable 250ml Bottle

IS1174 Oxygen Connecting Tube (1.8m)

SKF26-02 5 Year Service Pack

OVERVIEW

Oxylitre's range of Compact Medical Pipeline Flowmeters are the Flowmeter of choice for delivering flexible solutions. Our Compact Medical Pipeline Flowmeter Range are simple to use, reliable and extremely functional, our Compact Medical Pipeline Flowmeters offer high performance with durability engineered in to the robust design.

Oxylitre's range of Compact Medical Pipeline Flowmeters are dependable, extremely easy to maintain with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

