

COMPACT MEDICAL PIPELINE FLOWMETER

Compact Twin 'Y' Configuration Flowmeters Oxygen 0-15lpm Twin "Y" Configuration Flowmeters Oxygen 0-15Lpm Complete With Humidifier Bottles

PRODUCT CODE

F2602H



✓ 4 bar

✓ 0-15Lpm

SPECIFICATION

- OXYGEN
- Moulded BS Colour Coded Body
- NOT MRI Conditional
- "Y" Configuration Twin Mounting Block to British Standard Probe
- 9/16 UNF Fitting
- 4 bar nominal pressure
- BS Colour Fine adjustment Control
- Humidifier Bottle (Autoclavable to 121°C) complete with Tubing Nipple
- 0-15Lpm (O₂)

ACCESSORIES

F2602H

IS1174 Oxygen Connecting Tube (1.8m)

H110 Humidifier (Bubble Type) Autoclavable 250ml Bottle MRI Safe

SP46 Tubing Nipple female 9/16 UNF (White for O₂)

IS1106 Adult Face Mask O₂

OVERVIEW

Oxylitre's range of Compact Medical Pipeline Flowmeters are the Flowmeter of choice for delivering flexible solutions. Our Compact Medical Pipeline Flowmeter Range are simple to use, reliable and extremely functional, our Compact Medical Pipeline Flowmeters offer high performance with durability engineered in to the robust design.

Oxylitre's range of Compact Medical Pipeline Flowmeters are dependable, extremely easy to maintain with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom
+44 (0)161 872 6322 · sales@oxylitre.co.uk · service@oxylitre.co.uk · [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

