

COMPACT MEDICAL PIPELINE FLOWMETER

Compact Twin Flowmeter Medical Air 4 Bar 0-15Lpm Universal Rail Mounted

PRODUCT CODE

F2632A



✓ 4 Bar

✓ 0-15Lpm

SPECIFICATION

- Medical Air 4 Bar (MA-4)
- Moulded BS Colour Coded Body
- 9/16 UNF Fitting
- 4 bar nominal pressure
- BS Colour Fine adjustment Control
- "Y" Configuration Twin Mounting
- Tubing Nipple
- 0-15Lpm (MA-4)
- Complete with Universal Rail Clamp
- 1.5 metre BS Coloured Hose and Probe

ACCESSORIES

F2632A

H110 Humidifier (Bubble Type) Autoclavable 250ml Bottle

IS1174 Oxygen Connecting Tube (1.8m)

OVERVIEW

Oxylitre's range of Compact Medical Pipeline Flowmeters are the Flowmeter of choice for delivering flexible solutions. Our Compact Medical Pipeline Flowmeter Range are simple to use, reliable and extremely functional, our Compact Medical Pipeline Flowmeters offer high performance with durability engineered in to the robust design.

Oxylitre's range of Compact Medical Pipeline Flowmeters are dependable, extremely easy to maintain with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 · sales@oxylitre.co.uk · service@oxylitre.co.uk · oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29

MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

