

OXYGEN DEMAND UNIT

POD Oxygen Series Complete System 1.5m PVC Hose

PRODUCT CODE

POD2002



SPECIFICATION

- POD Oxygen Series Complete System
- Demand Head complete with 1.5m Hose and BS Probe
- Antimicrobial Infused Demand Head
- 360° fully rotational integral handle
- Exhaust Vent away from user
- Uses In-line filters
- Compact Regulator Self Sealing Valve Pin-Index
- Blue Weatherproof Antimicrobial Carry Case
- 1 x Pin-Index Key & Cord
- Consumable starter kit consisting of:
 - 1 x Adult Face Mask
 - 1 x Child Face Mask
 - 1 x Bacterial Filter Disk
 - 1 x Pack of 5 Mouth pieces

ACCESSORIES

POD2002

POK050 Oxygen Demand System Carry Case

PRU010 Bacterial Fliter Disks (Box of 50)

332555 Mouthpieces (Packet of 5)

SKP23 5 Year Service Pack

OVERVIEW

Oxylitre's range of Oxygen Demand Units are the oxygen demand valve of choice. Delivering numerous and significant improvements whilst driving down the cost of a high quality dependable oxygen demand valve. The Demand Valve is available in a full range of sets, designed to meet your needs from field use to hospital, Oxylitre offer a oxygen therapy solution for you.

Oxylitre Oxygen Demand Units are user friendly, intuitive to hold with a 360 degree fully rotational integral handle. The Demand Unit has been comprehensively field tested made from high impact plastic, the units streamline compact design is both durable and lightweight. Our Oxygen Demand Valve is trouble free and easy to maintain with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 · sales@oxylitre.co.uk · service@oxylitre.co.uk · oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

