

ENTONOX® PAIN RELIEF SYSTEM

Entonox® Pain Relief System Complete 3m PVC Hose

PRODUCT CODE

PRU2002L



SPECIFICATION

- Entonox® Pain Relief System Complete Mk2
- Demand Head complete with 3m Hose and BS Probe
- Antimicrobial Infused Demand Head
- 360° fully rotational integral handle
- Exhaust Vent away from user
- Uses In-line filters
- Compact Regulator Self Sealing Valve Pin-Index
- Blue Weatherproof Antimicrobial Carry Case
- 1 x Pin-Index Key & Cord
- Consumable starter kit consisting of:
 - 1 x Adult Face Mask
 - 1 x Child Face Mask
 - 1 x Bacterial Filter Disk
 - 1 x Pack of 5 Mouth pieces

ACCESSORIES

PRU2002L

PRU010 Bacterial Filter Disks (Box of 50)

332555 Mouthpieces (Packet of 5)

SKP23 5 Year Service Pack

OVERVIEW

Oxylitre' range of Pain Relief Units (the PRU), are the demand valve of choice. Delivering numerous and significant improvements whilst driving down the cost of a high quality dependable emergency demand valve. The Demand Valve is available in a full range of sets, designed to meet your needs from field use to maternity. Oxylitre offer a Pain Relief solution for you.

The PRU' is user friendly, intuitive to hold with a 360 degree fully rotational integral handle. The PRU has been comprehensively field tested made from high impact plastic, the units streamline compact design is both durable and lightweight. The PRU' is trouble free and easy to maintain with a long service life, with a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

