

## PORTABLE SUCTION PUMP

# "petite elite" Portable Suction Pump High Suction VaxSax Jar & Liner (Portable)

PRODUCT CODE

**PSP004**



✓ 500mmHg

✓ 25Lpm

✓ 110volt

✓ 12volt

## SPECIFICATION

- High Vacuum High Flow
- Single 1000ml Disposable VACSAX Jar & Liner System
- Robust High Impact Resistant Moulded Enclosure
- Encapsulated Bacterial/Hydrophobic Filter Disk
- Luminous Touch Panel Control
- Instant Maximum Vacuum Override Button
- Battery Level / Charging Indicator
- 80 to 500mmHg.
- Flowrate >25Lpm
- 110volt – 240volt AC
- 12volt – 28volt DC
- 12 volt DC Rechargeable Battery, 4.3mAh
- Run time >1hr 20min
- Weight 4.9Kg
- Approximate Size 310mm x 305mm x 170mm
- Peak Sound Output at MAX Vacuum from 1m: 65 decibels

## WEIGHTS AND DIMENSIONS

|        |         |
|--------|---------|
| Weight | 4.43 kg |
| Height | 320 mm  |
| Width  | 300 mm  |
| Depth  | 170 mm  |

## ACCESSORIES

PSP004

**PSP600-50** Bacterial/Hydrophobic Filter Disks (box of 50)

**SKPSP5** Major Service Pack

## OVERVIEW

Oxylitre's range of Portable Suction Pumps (PSP's) deliver a high quality flexible portable suction solution for Healthcare and Emergency Services. Our contemporary design Portable Suction is user friendly, comprehensively field tested to the latest standards. Made from high impact plastic, the units streamline compact design is both durable and lightweight. The Portable Suction Pumps (PSP's) are intuitive to use incorporating a luminous LED illuminated touch control panel and an easy access battery compartment. These units are reliable and extremely functional our range of Portable Suction Units offer high performance with operational reliability.

Oxylitre's range of Portable Suction Units are dependable, with minimum maintenance required. Our Range of Portable Suction Units come with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom  
+44 (0)161 872 6322 · sales@oxylitre.co.uk · service@oxylitre.co.uk · [oxylitre.co.uk](https://oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29  
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

