

STANDARD MEDICAL REGULATOR

Regulator Twin Gauge Medical Air 0-15lpm

PRODUCT CODE

R15005A



✓ 138 bar

✓ 4 bar

✓ 0-15Lpm

SPECIFICATION

- Medical Air
- Brass Chrome Plated Body
- NOT MRI Conditional
- ISO 407 – Pin-Index Cylinder Fitting
- 138 bar maximum inlet pressure
- Pre-Set 4 bar output pressure
- High Pressure Cylinder Content Gauge
- Complete with Gauge Flowmeter:
- Tubing Nipple
- 0-15Lpm (MA-4)

WEIGHTS AND DIMENSIONS

Weight	1.52 kg
Height	179 mm
Width	139 mm
Depth	191 mm

ACCESSORIES

R15005A

IS1174 Oxygen Connecting Tube (1.8m)

BS110 Bullnose O' Ring Standard Type (each)

OVERVIEW

Oxylitre's range of Standard Medical Regulators, have been designed with the end user in mind. Our Standard Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Regulators offer high performance with durability engineered in to the robust design.

Oxylitre's range of Standard Medical Regulators are dependable, easy to maintain with a long service life with a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom
+44 (0)161 872 6322 · sales@oxylitre.co.uk · service@oxylitre.co.uk · [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

