

STANDARD MEDICAL REGULATOR

Regulator & Flowmeter Carbon Dioxide 0-3lpm

PRODUCT CODE

R1600-C35



✓ 50 bar

✓ 4 bar

✓ 0-3Lpm

SPECIFICATION

- Carbon Dioxide
- Brass Chrome Plated Body
- NOT MRI Conditional
- ISO 407 – Pin-Index Cylinder Fitting
- 50 bar inlet pressure
- Pre-Set 4 bar output pressure
- Cylinder content Gauge
- Complete with Flowmeter:
- Flowmeter Brass Chrome Plated 9/16 UNF
- 3/8" Nut & Liner Connection
- Tubing Nipple
- 0-3Lpm (CO₂)

WEIGHTS AND DIMENSIONS

Weight	1.49 kg
Height	243 mm
Width	181 mm
Depth	194 mm

ACCESSORIES

R1600-C35

SP46 Tubing Nipple female 9/16 UNF (White for O₂)

IS1174 Oxygen Connecting Tube (1.8m)

BS110 Bullnose O' Ring Standard Type (each)

OVERVIEW

Oxylitre's range of Standard Medical Regulators, have been designed with the end user in mind. Our Standard Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Regulators offer high performance with durability engineered in to the robust design.

Oxylitre's Standard Medical Regulators are dependable and easy to maintain, with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom
+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

