

STANDARD MEDICAL REGULATOR

Regulator Pre-Set Nitrous Oxide Left-Handed 11/16" Bullnose Fitting

PRODUCT CODE

R1605LN-BN



✓ 44 bar

✓ 4 bar

SPECIFICATION

- Nitrous Oxide
- Brass Chrome Plated Body
- NOT MRI Conditional
- BS 341 No. 13 11/16" Male Cylinder Fitting
- 44 bar maximum inlet pressure
- Pre-Set 4 bar output pressure
- 1/8 bsp Left Handed Male Non Return Connection

WEIGHTS AND DIMENSIONS

Weight	0.8 kg
Height	88 mm
Width	89 mm
Depth	68 mm

ACCESSORIES

R1605LN-BN

SP46 Tubing Nipple female 9/16 UNF (White for O2)

IS1174 Oxygen Connecting Tube (1.8m)

BS110 Bullnose O' Ring Standard Type (each)

OVERVIEW

Oxylitre's range of Standard Medical Regulators, have been designed with the end user in mind. Our Standard Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Regulators offer high performance with durability engineered in to the robust design.

Oxylitre's range of Standard Medical Regulators are dependable, easy to maintain with a long service life with a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29

MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

