

STANDARD MEDICAL REGULATOR

Regulator 11/16" Nut & Tubing Outlet Nitrous Oxide

PRODUCT CODE

R16105N



✓ 44 bar

✓ 4 bar

SPECIFICATION

- Nitrous Oxide
- Brass Chrome Plated Body
- NOT MRI Conditional
- Pin-Index Cylinder Fitting (Compatible with UK Size D & E Nitrous Oxide Cylinders)
- 44 bar maximum inlet pressure
- Pre-Set 4 bar output pressure
- Cylinder content Gauge
- N₂O 11/16" Nut & Tubing Outlet Connection

WEIGHTS AND DIMENSIONS

Weight	1.29 kg
Height	189 mm
Width	132 mm
Depth	177 mm

ACCESSORIES

R16105N

SK60 5 Year Service Kit

BS110 Bullnose O' Ring Standard Type (each)

OVERVIEW

Oxylitre's range of Standard Medical Regulators, have been designed with the end user in mind. Our Standard Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Regulators offer high performance with durability engineered in to the robust design.

Oxylitre's range of Standard Medical Regulators are dependable, easy to maintain with a long service life with a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

