

STANDARD MEDICAL REGULATOR

Regulator & Self Sealing Valve Medical Air (MA-7)

PRODUCT CODE

R1611



✓ 7 Bar

✓ 138 bar

SPECIFICATION

- Medical Air 7 Bar (MA-7)
- Brass Chrome Plated Body
- NOT MRI Conditional
- High Pressure Cylinder Content Gauge
- Bullnose Cylinder Fitting (Compatible with UK Size F & G Medical Air Cylinders)
- 138 bar maximum inlet pressure
- Pre-Set 7 bar output pressure
- British Standard Gas Indexed Coloured Self Sealing Valve

WEIGHTS AND DIMENSIONS

Weight	1.22 kg
Height	171 mm
Width	124 mm
Depth	74 mm

ACCESSORIES

R1611

SK60 5 Year Service Kit

BS110 Bullnose O' Ring Standard Type (each)

OVERVIEW

Oxylitre's range of Standard Medical Regulators, have been designed with the end user in mind. Our Standard Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Regulators offer high performance with durability engineered in to the robust design.

Oxylitre's range of Standard Medical Regulators are dependable, easy to maintain with a long service life with a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 · sales@oxylitre.co.uk · service@oxylitre.co.uk · [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

