

## STANDARD MEDICAL REGULATOR

# Regulator & Self Sealing Valve Medical Air (MA-7)

PRODUCT CODE

**R16115**



✓ 7 Bar

✓ 138 bar

## SPECIFICATION

- Medical Air 7 Bar (MA-7)
- Brass Chrome Plated Body
- NOT MRI Conditional
- High Pressure Cylinder Content Gauge
- ISO 407 - Pin-Index Cylinder Fitting
- 138 bar maximum inlet pressure
- Pre-Set 7 bar output pressure
- British Standard Gas Indexed Coloured Self Sealing Valve

## WEIGHTS AND DIMENSIONS

|        |         |
|--------|---------|
| Weight | 1.43 kg |
| Height | 191 mm  |
| Width  | 124 mm  |
| Depth  | 190 mm  |

## ACCESSORIES

R16115

**SK60** 5 Year Service Kit

**BS110** Bullnose O' Ring Standard Type (each)

## OVERVIEW

Oxylitre's range of Standard Medical Regulators, have been designed with the end user in mind. Our Standard Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Regulators offer high performance with durability engineered in to the robust design.

Oxylitre's range of Standard Medical Regulators are dependable, easy to maintain with a long service life with a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29  
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

