


STANDARD MEDICAL REGULATOR

Regulator High Pressure Medical Air (MA-7)

PRODUCT CODE

R17005



✓ 7 Bar

✓ 138 bar

SPECIFICATION

- Medical Air 7 Bar (MA-7)
- Brass Chrome Plated Body
- NOT MRI Conditional
- ISO 407 – Pin-Index Cylinder Fitting
- 138 bar maximum inlet pressure
- Pre-Set 7 bar output pressure
- High Pressure Cylinder Content Gauge
- Low Pressure Gauge 0 to 150 PSI
- Adjustable Output Pressure
- British Standard Gas Indexed Coloured Self Sealing Valve

WEIGHTS AND DIMENSIONS

Weight	1.77 kg
Height	163 mm
Width	156 mm
Depth	228 mm

ACCESSORIES

R17005

AH103-3 Extention Hose Assembly BS Probe to SSV

BS110 Bullnose O' Ring Standard Type (each)

OVERVIEW

Oxylitre's range of Standard Medical Regulators, have been designed with the end user in mind. Our Standard Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Regulators offer high performance with durability engineered in to the robust design.

Oxylitre's range of Standard Medical Regulators are dependable, easy to maintain with a long service life with a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 · sales@oxylitre.co.uk · service@oxylitre.co.uk · [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

