

## STANDARD MEDICAL REGULATOR

# Regulator High Pressure Medical Air (MA-7) 3/8bsp Outlet

PRODUCT CODE

**R1700A**



✓ 7 Bar

✓ 138 bar

## SPECIFICATION

- Medical Air 7 Bar (MA-7)
- Brass Chrome Plated Body
- NOT MRI Conditional
- Bullnose Cylinder Fitting (Compatible with UK Size F & G Cylinders – BS 341-3)
- 138 bar maximum inlet pressure
- Pre-Set 7 bar output pressure
- High Pressure Cylinder Content Gauge
- Low Pressure Gauge 0 to 150 PSI
- Adjustable Output Pressure
- 3/8bsp Male Outlet

## WEIGHTS AND DIMENSIONS

|        |         |
|--------|---------|
| Weight | 1.49 kg |
| Height | 154 mm  |
| Width  | 138 mm  |
| Depth  | 135 mm  |

## ACCESSORIES

R1700A

**AH103-3** Extention Hose Assembly BS Probe to SSV

**BS110** Bullnose O' Ring Standard Type (each)

## OVERVIEW

Oxylitre's range of Standard Medical Regulators, have been designed with the end user in mind. Our Standard Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Regulators offer high performance with durability engineered in to the robust design.

Oxylitre's range of Standard Medical Regulators are dependable, easy to maintain with a long service life with a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29

MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

