

COMPACT MEDICAL REGULATOR

Compact Regulator & Brass Flowmeter Oxygen 0-15lpm c/w Humidifier Bottle

PRODUCT CODE

R24005H



✓ 200 bar

✓ 4 bar

✓ 0-15Lpm

SPECIFICATION

- OXYGEN
- Stainless Steel Body
- ISO 407 – Pin-Index Cylinder Fitting
- 200 bar maximum inlet pressure
- Pre-Set 4 bar output pressure
- High Pressure Luminous Cylinder Content Gauge
- Complete with Flowmeter:
- Flowmeter Brass Chrome Plated 9/16 UNF
- Humidifier Bottle (Autoclavable to 121°C) complete with Tubing Nipple
- 0-15Lpm (O₂)

WEIGHTS AND DIMENSIONS

Weight	0.94 kg
Height	320 mm
Width	150 mm
Depth	75 mm

ACCESSORIES

R24005H

IS1174 Oxygen Connecting Tube (1.8m)

SK40 Service Kit for R14 & R24 Regulator Range (Bullnose & Pin-Index)

H110 Humidifier (Bubble Type) Autoclavable 250ml Bottle MRI Safe

BS110 Bullnose O' Ring Standard Type (each)

OVERVIEW

Oxylitre's range of Compact Medical Regulators have been designed with the end user in mind. Rated up to 200 bar input pressure, our Compact Range of regulators are simple to use, reliable and extremely functional, offering high performance with durability engineered in to the robust design.

Oxylitre's Compact Medical Regulators are dependable and easy to maintain, with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

