

COMPACT MEDICAL REGULATOR

Compact Regulator & Brass Flowmeter Air 0-15lpm

PRODUCT CODE

R2400A



✓ 4 bar

✓ 200 bar

✓ 0-15Lpm

SPECIFICATION

- Medical Air 4 bar (MA-4)
- Stainless Steel Body
- Bullnose Cylinder Fitting (Compatible with UK Size F & G Medical Air Cylinders)
- 200 bar maximum inlet pressure
- Pre-Set 4 bar output pressure
- High Pressure Luminous Cylinder Content Gauge
- Complete with Flowmeter:
- Flowmeter Brass Chrome Plated 9/16 UNF
- 3/8" Nut & Liner Connection
- Air Safety Shield & Tubing Nipple
- 0-15Lpm (MA-4)

WEIGHTS AND DIMENSIONS

Weight	0.84 kg
Height	216 mm
Width	148 mm
Depth	62 mm

ACCESSORIES

R2400A

SP46 Tubing Nipple female 9/16 UNF (White for O2)

IS1174 Oxygen Connecting Tube (1.8m)

BS110 Bullnose O' Ring Standard Type (each)

SK40 5 Year Service Pack

OVERVIEW

Oxylitre's range of Compact Medical Regulators have been designed with the end user in mind. Rated up to 200 bar input pressure, our Compact Range of regulators are simple to use, reliable and extremely functional, offering high performance with durability engineered in to the robust design.

Oxylitre's Compact Medical Regulators are dependable and easy to maintain, with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

