

COMPACT MEDICAL REGULATOR

Compact Regulator & Moulded Flowmeter Oxygen 0-15lpm c/w Humidifier Bottle

PRODUCT CODE

R2400MH



✓ 200 bar

✓ 4 bar

✓ 0-15Lpm

SPECIFICATION

- OXYGEN
- Stainless Steel Body
- Bullnose Cylinder Fitting (Compatible with UK Size F & G Cylinders – BS 341–3)
- 200 bar maximum inlet pressure
- Pre-Set 4 bar output pressure
- High Pressure Luminous Cylinder Content Gauge
- Complete with Flowmeter:
- Moulded Flowmeter 9/16 UNF
- Humidifier Bottle (Autoclavable to 121°C) complete with Tubing Nipple
- 0-15Lpm (O₂)

WEIGHTS AND DIMENSIONS

Height	330 mm
Width	148 mm
Depth	62 mm

ACCESSORIES

R2400MH

SP46 Tubing Nipple female 9/16 UNF (White for O₂)

IS1174 Oxygen Connecting Tube (1.8m)

BS110 Bullnose O' Ring Standard Type (each)

SK40 5 Year Service Pack

OVERVIEW

Oxylitre's range of Compact Medical Regulators have been designed with the end user in mind. Rated up to 200 bar input pressure, our Compact Range of regulators are simple to use, reliable and extremely functional, offering high performance with durability engineered in to the robust design.

Oxylitre's Compact Medical Regulators are dependable and easy to maintain, with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom
+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

