

COMPACT SELECT-A-FLOW REGULATOR

Compact Select-A-Flow Regulator 3/8" BSP Outlet Pin-Index Oxygen 0-15lpm

PRODUCT CODE

R24035-10



✓ 200 bar

✓ 4 bar

✓ 0-15Lpm

SPECIFICATION

- OXYGEN
- Stainless Steel Body
- ISO 407 – Pin-Index Cylinder Fitting
- 200 bar maximum inlet pressure
- Pre-Set 4 bar output pressure
- High Pressure Luminous Cylinder Content Gauge
- Complete with Flowmeter:
- Select-A-Flow Flowmeter (11 Settings)(Other flow rates available)
- 3/8 BSP Male Outlet
- 0-15Lpm (O₂)

WEIGHTS AND DIMENSIONS

Weight	0.76 kg
Height	193 mm
Width	65 mm
Depth	72 mm

ACCESSORIES

R24035-10

IS1174 Oxygen Connecting Tube (1.8m)

BS110 Bullnose O' Ring Standard Type (each)

SK40 5 Year Service Pack

OVERVIEW

Oxylitre's range of Compact Medical Regulators have been designed with the end user in mind. Rated up to 200 bar input pressure, our Compact Range of regulators are simple to use, reliable and extremely functional, offering high performance with durability engineered in to the robust design.

Oxylitre's range of Compact Medical Regulators are dependable, easy to maintain with a long service life with a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

