

COMPACT MEDICAL REGULATOR

Compact Select-A-Flow Regulator Bullnose Oxygen 0-15lpm & Self Sealing Valve

PRODUCT CODE

R2405



✓ 200 bar

✓ 4 bar

✓ 0-15Lpm

SPECIFICATION

- OXYGEN
- Stainless Steel Body
- Bullnose Cylinder Fitting (Compatible with UK Size F & G Cylinders – BS 341–3)
- 200 bar maximum inlet pressure
- Pre-Set 4 bar output pressure
- High Pressure Luminous Cylinder Content Gauge
- Complete with Flowmeter:
- Select-A-Flow Flowmeter (11 Settings)(Other flow rates available)
- Tubing Nipple
- 0-15Lpm (O₂)
- British Standard Gas Indexed Coloured Self Sealing Valve

WEIGHTS AND DIMENSIONS

Weight	0.90 kg
Height	151 mm
Width	93 mm
Depth	72 mm

ACCESSORIES

R2405

IS1174 Oxygen Connecting Tube (1.8m)

OX010 Pin-Index "BODOK" Seal (each)

SK40 5 Year Service Pack

OVERVIEW

Oxylitre's range of Compact Medical Regulators have been designed with the end user in mind. Rated up to 200 bar input pressure, our Compact Range of regulators are simple to use, reliable and extremely functional, offering high performance with durability engineered in to the robust design.

Oxylitre's range of Compact Medical Regulators are dependable, easy to maintain with a long service life with a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom
+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

