


COMPACT MEDICAL REGULATOR

Compact Regulator 3/8bsp Male Outlet Oxygen

PRODUCT CODE

R24105



✓ 200 bar

✓ 4 bar

SPECIFICATION

- OXYGEN
- Stainless Steel Body
- ISO 407 – Pin-Index Cylinder Fitting
- 200 bar maximum inlet pressure
- Pre-Set 4 bar output pressure
- High Pressure Luminous Cylinder Content Gauge
- 3/8" BSP Male Outlet - Nut & Liner Connection

WEIGHTS AND DIMENSIONS

Weight	0.62 kg
Height	143 mm
Width	61 mm
Depth	62 mm

ACCESSORIES

R24105

SP46 Tubing Nipple female 9/16 UNF (White for O2)

IS1174 Oxygen Connecting Tube (1.8m)

BS110 Bullnose O' Ring Standard Type (each)

SK40 5 Year Service Pack

OVERVIEW

Oxylitre's range of Compact Medical Regulators have been designed with the end user in mind. Rated up to 200 bar input pressure, our Compact Range of regulators are simple to use, reliable and extremely functional, offering high performance with durability engineered in to the robust design.

Oxylitre's range of Compact Medical Regulators are dependable, easy to maintain with a long service life with a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

