

ON/OFF MEDICAL PIPELINE SUCTION

"elite" Twin Pipeline Suction Controllers High/Thoracic Suction Direct British Standard Probe

PRODUCT CODE

S723-HT



SPECIFICATION

- High Vacuum High Flow
- Thoracic Vacuum
- Robust Moulded Construction
- British Standard Probe
- Integral Gauge colour coded 0-760mmHg: DUAL Scale: kPa & mmHg
- Integral Gauge colour coded 0-200mmHg: DUAL Scale: kPa & mmHg
- Positive Pressure Relief Valve
- On/Off Flag Style Lever
- Fine adjustment On Off Increase/Decrease Control
- Encapsulated Bacterial/Hydrophobic Filter Unit

WEIGHTS AND DIMENSIONS

Weight	1.14 kg
Height	175 mm
Width	220 mm
Depth	190 mm

ACCESSORIES

S723-HT

- S750** Hydrophobic/Bacterial Filters (Box of 30)
- S7500A** 2 Litre Suction Jar Assembly "V" Mounted
- S7500B** 2 Litre Suction Jar Assembly Rail Mounted
- SK17** 5 Year Service Pack

OVERVIEW

Oxylitre's range of Medical Pipeline Units deliver flexible range of suction solutions. Our Medical Pipeline Suction Range have a long and exemplary history, simple to use, reliable and extremely functional. The Medical Pipeline Suction Range offers high performance with operational reliability engineered in to the robust design.

Oxylitre's range of Medical Pipeline Suction are dependable, with minimum maintenance required. Our Range of Portable Suction Units come with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom
+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

