

ON/OFF MEDICAL PIPELINE SUCTION

"elite" Pipeline Suction Controller Thoracic Suction Remote Universal Rail Clamp

PRODUCT CODE

S725R



SPECIFICATION

- Universal Rail Clamp
- 1.5 metre BS Coloured Hose and BS Probe connection
- Robust Moulded Construction
- Integral Gauge colour coded 0-760mmHg: DUAL Scale: kPa & mmHg
- Positive Pressure Relief Valve
- On/Off Flag Style Lever
- Fine adjustment On Off Increase/Decrease Control
- Encapsulated Bacterial/Hydrophobic Filter Unit

WEIGHTS AND DIMENSIONS

Weight	0.93 kg
Height	175 mm
Width	92 mm
Depth	160 mm

ACCESSORIES

S725R

- S750** Hydrophobic/Bacterial Filters (Box of 30)
- S7500A** 2 Litre Suction Jar Assembly "V" Mounted
- S7500B** 2 Litre Suction Jar Assembly Rail Mounted
- SK17** 5 Year Service Pack

OVERVIEW

Oxylitre's range of Medical Pipeline Units deliver flexible range of suction solutions. Our Medical Pipeline Suction Range have a long and exemplary history, simple to use, reliable and extremely functional. The Medical Pipeline Suction Range offers high performance with operational reliability engineered in to the robust design.

Oxylitre's range of Medical Pipeline Suction are dependable, with minimum maintenance required. Our Range of Portable Suction Units come with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

