

SUCTION JAR ASSEMBLY

2 Litre Suction Jar Assembly Rail Mounted

PRODUCT CODE

S7500B



SPECIFICATION

- 2 Litre Capacity Jar
- Autoclavable to 134°C
- Overfill Float Protection
- Stainless Steel Cradle
- Universal Rail Clamp

WEIGHTS AND DIMENSIONS

Weight	0.7 kg
Height	300 mm
Width	200 mm
Depth	150 mm

ACCESSORIES

S7500B

S750 Hydrophobic/Bacterial Filters (Box of 30)

S7500B 2 Litre Suction Jar Assembly Rail Mounted

OVERVIEW

Oxylitre's range of Medical Pipeline Units deliver flexible range of suction solutions. Our Medical Pipeline Suction Range have a long and exemplary history, simple to use, reliable and extremely functional. The Medical Pipeline Suction Range offers high performance with operational reliability engineered in to the robust design.

Oxylitre's range of Standard Medical Pipeline Suction are dependable, with a minimum maintain required. Our Range of Medical Pipeline Suction Units come with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

