

INJECTOR SUCTION SYSTEMS

Patient Oxygen & Suction Transfer System (POST)

PRODUCT CODE

S8000



✓ 600mmHg

✓ 4 bar

SPECIFICATION

	Colour coordinated high vacuum injector suction controller
	Variable vacuum 0 to >600mmHg
	Colour coordinated "click style" Select-a-Flow flowmeter
	Vertical or horizontal mounting positions
	Rail clamp
	1 metre hose and probe
	4 bar positive pressure input
Select-a-Flow Settings	0 to 15 Lpm: 0, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8, 10, 15

ACCESSORIES

S8000

- SK8000** Service Kit S8000 P.O.S.T Unit
- S8701A** Cradle & Rail Mount Assembly
- VAC100** Suction Tubing Clear A/S (30m Roll)
- IV008** Injector Silencer Wadding (2 per unit) each

OVERVIEW

Oxylitre's range of Injector Suction Systems deliver practical, efficient suction solutions for healthcare services. Our Patient Transfer Suction Systems have a long and exemplary history: simple to use, reliable and extremely functional. They are compact and extremely durable, designed for hard working environments and delivering a repeatable high performance time after time.

Oxylitre's range of Patient Transfer Injector Suction Systems are dependable and require minimum maintenance. Our range of Patient Transfer Injector Suction Systems come with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

