

PATIENT TRANSFER INJECTOR SUCTION SYSTEM

Patient Transfer Injector Suction System Basic (Flowmeter & 3/8bsp connector)

PRODUCT CODE

S801-1



SPECIFICATION

- High Vacuum Injector Suction Controller 0 to >500mmHg
- Flowmeter 0-15Lpm
- Back Bar Mounting Assembly
- 3/8 BSP Nut and Liner Connection (To fit Regulator)

ACCESSORIES

S801-1

VAC201 Yellow Suction Tubing 1m

S7500A 2 Litre Suction Jar Assembly "V" Mounted

S7500B 2 Litre Suction Jar Assembly Rail Mounted

OVERVIEW

Oxylitre's range of Injector Suction Systems deliver practical efficient suction solutions for Healthcare Services. Oxylitre's range of Patient Transfer Suction Systems have a long and exemplary history, simple to use, reliable and extremely functional. Our Patient Transfer Suction Systems are compact and extremely durable designed for hard working environments delivering a repeatable high performance time after time.

Oxylitre's range of Patient Transfer Injector Suction Systems are dependable and require minimum maintenance. Our Range of Patient Transfer Injector Suction Systems come with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

