

COMPACT POST SYSTEM

1 Litre Suction Jar Assembly "V" Mounted

PRODUCT CODE

S8702



SPECIFICATION

- 1000ml Autoclavable to 121°C Suction Jar Assembly
- Screw on Lid complete with Tube connections
- Overflow Float Assembly, Autoclavable
- Stainless Steel Cradle
- Jar Retaining Clips "V" Mounted Assembly

ACCESSORIES

S8702

VAC201 Suction Tubing Yellow (per 1m)

OVERVIEW

Oxylitre's range of Injector Suction Systems deliver practical efficient suction solutions for Healthcare Services. Oxylitre's range of Patient Transfer Suction Systems have a long and exemplary history, simple to use, reliable and extremely functional. Our Patient Transfer Suction Systems are compact and extremely durable designed for hard working environments delivering a repeatable high performance time after time.

Oxylitre's range of Patient Transfer Injector Suction Systems are dependable and require minimum maintenance. Our Range of Patient Transfer Injector Suction Systems come with a long service life and a 7 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

