

## PIPELINE AND HOSE ACCESSORIES

# British Standard High Pressure Hose

PRODUCT CODES  
**HOSE1 to HOSE7**



✓ 12.7mm

✓ 6.3mm

## SPECIFICATION

- ESD Anti-Static Inner Core
- Antimicrobial Outer Core
- Non-Phthalate PVC - DEHP Free
- Reach Compliant - SVHC Free
- Hose O/D 12.7mm
- Hose I/D 6.3mm
- Standard roll 30m

## RANGE

**HOSE1** White - Oxygen

**HOSE2** Blue - Nitrous Oxide

**HOSE3** Black - Medical Air

**HOSE4** Yellow - Vacuum

**HOSE6** White/Blue - Entonox®

**HOSE7** Grey - Carbon Dioxide

## OVERVIEW

Oxylitre Holdings Limited offer the most comprehensive range of Pipeline and Hose Accessories. Our Pipeline and Hose Accessories cover all Medical Gases suitable. Our connections are suitable for use with all flexible medical hose. Our extensive range of gas specific pipeline connections and accessories cover both British and European Standard connections. A full and comprehensive look at our range can be found in our product index.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom  
+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29  
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

