

PORTABLE SUCTION PUMPS

Portable Suction Pump Docking Stations

PRODUCT CODES
6 codes



✓ 1.5Kg

✓ 70mm

SPECIFICATION

- Lightweight (1.5Kg approx)
- Automatic Unit Charging (when powered)
- Tested Beyond Vehicle Transit Requirements (20g)
- Approx Size 310 x 305 x 70mm

RANGE

- PSP701A** 12v DC with 1.5m Lead
- PSP701B** 24-28v DC with 1.5m Lead
- PSP701C** 240v AC Lead c/w UK 3 Pin Plug, 1.5m Lead
- PSP701D** 240v AC Lead c/w EU 2 Pin Plug, 1.5m Lead
- CSP701** Docking Station Only (without power supply)
- PSP701M** Docking Station Mounting Bush Kit (3 Hole)

ACCESSORIES

6 CODES

- PSP600-50** Bacterial/Hydrophobic Filter Disks (Box of 50)
- PSP500** 1000ml STANDARD Jar Assembly complete
- PSP500-5-1** Suction Catheter Connector (each)

OVERVIEW

Oxylitre's range of Portable Suction Pumps (PSP's) deliver a high quality flexible portable suction solution for Healthcare and Emergency Services. Our contemporary design Portable Suction is user friendly, comprehensively field tested to the latest standards. Made from high impact plastic, the units streamline compact design is both durable and lightweight. The Portable Suction Pumps (PSP's) are intuitive to use incorporating a luminous LED illuminated touch control panel and an easy access battery compartment. These units are reliable and extremely functional our range of Portable Suction Units offer high performance with operational reliability.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom
+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

