

PIPELINE AND HOSE ACCESSORIES

Pipeline Connectors/Fittings

PRODUCT CODES

9 codes



SPECIFICATION

- See Product Range

RANGE

SP70 Nut 3/8bsp (Right Hand)

SP70A Tubing Outlet for 3/8bsp Nut

SP75 Nut and Tubing Tail 1/4bsp

SB30 Nut 3/8bsp (Left Hand)

SP92A O₂ DISS Nut (9/16unf)

SP90A O₂ DISS Liner (Plain) x 6.5mm O/D Tubing Tail

SP90 O₂ DISS Liner (Grooved) x 6.5mm O/D Tubing Tail

6363 Connector 1/4bsp x 6.5mm O/D Tubing Tail

BS010 British Standard 6.07mm x 1.78mm O' Ring

OVERVIEW

Oxylitre Holdings Limited offer the most comprehensive range of Pipeline and Hose Accessories. Our Pipeline and Hose Accessories cover all Medical Gases suitable. Our connections are suitable for use with all flexible medical hose. Our extensive range of gas specific pipeline connections and accessories cover both British and European Standard connections. A full and comprehensive look at our range can be found in our product index.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

