

PIPELINE HOSE ASSEMBLY NIST FITTING TO SELF SEALING VALVE

Hose Assy Oxygen, 3m, NIST Nut to S/S Valve

PRODUCT CODES

AH201 to AH206



SPECIFICATION

- British Standard NIST Coupling
- British Standard Coloured Indexed Self Sealing Valve
- 3 Meter British Standard Colour Coded Hose
- Anti Static and Phthalate Free

RANGE

- AH201** Hose Assy Oxygen, 3m, NIST Nut to S/S Valve
- AH202** Hose Assy N2O, 3m, NIST Nut to S/S Valve
- AH203** Hose Assy Entonox, 3m, NIST Nut to S/S Valve
- AH204** Hose Assy MA-4, 3m, NIST Nut to S/S Valve
- AH205** Hose Assy MA-7, 3m, NIST Nut to S/S Valve
- AH206** Hose Assy VACUUM, 3m, NIST Nut to S/S Valve

WEIGHTS AND DIMENSIONS

Dimensions

L - 3m

OVERVIEW

Oxylitre Holdings Limited offer the most comprehensive range of Pipeline Hose Assemblies. Our Pipeline Hose Assemblies Range cover all Medical Gases from oxygen through to special gas mixtures. Made from flexible anti-static phalate free colour coordinated hose. Oxylitre also offer a full design service, to create the hose assembly you require, for more information on this service please contact sales. Our listed hose assemblies comprise of an extensive range of gas specific connections. A full and comprehensive look at our range can be found in our product index.

Each Pipeline Hose Assembly is produced to BS and European Standards and all are available in none listed lengths. Oxylitre's Hose Assembly Range comes with a minimum 5 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29

MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

