

EXTENSION HOSE ASSEMBLY BRITISH STANDARD PROBE TO SELF SEALING VALVE

Hose Assy Extension Oxygen, 3m, BS Probe to S/S Valve

PRODUCT CODES

AH101-3 to AH106-3



SPECIFICATION

- British Standard Probe
- British Standard Coloured Indexed Self Sealing Valve
- 3 Meter British Standard Colour Coded Hose
- Anti Static and Phthalate Free

RANGE

AH101-3 Hose Assy Extension Oxygen, 3m, BS Probe to S/S Valve

AH102-3 Hose Assy Extension N2O, 3m, BS Probe to S/S Valve

AH103-3 Hose Assy Extension MA-7, 3m, BS Probe to S/S Valve

AH104-3 Hose Assy Extension MA-4, 3m, BS Probe to S/S Valve

AH105-3 Hose Assy Extension VACUUM, 3m, BS Probe to S/S Valve

AH106-3 Hose Assy Extension Entonox, 3m, BS Probe to S/S Valve

WEIGHTS AND DIMENSIONS

Dimensions

L - 3m

OVERVIEW

Oxylitre Holdings Limited offer the most comprehensive range of Pipeline Hose Assemblies. Our Pipeline Hose Assemblies Range cover all Medical Gases from oxygen through to special gas mixtures. Made from flexible anti-static phalate free colour coordinated hose. Oxylitre also offer a full design service, to create the hose assembly you require, for more information on this service please contact sales. Our listed hose assemblies comprise of an extensive range of gas specific connections. A full and comprehensive look at our range can be found in our product index.

Each Pipeline Hose Assembly is produced to BS and European Standards and all are available in none listed lengths. Oxylitre's Hose Assembly Range comes with a minimum 5 year warranty.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

