

PIPELINE ACCESSORIES

Oxygen 4 Bar Input (42mm wide)

PRODUCT CODES

T1000R to T1005RW



SPECIFICATION

- Brass Chrome Plated Body
- NOT MRI Conditional
- BS Indexed Coloured Self Sealing Valve
- 1.5m British Standard Hose and Probe Assembly
- Universal Rail Clamp

RANGE

- T1000R** Oxygen 4 Bar Input (42mm wide)
- T1000RW** Oxygen 4 Bar Input (94mm wide)
- T1001R** Nitrous Oxide 4 Bar Input (42mm wide)
- T1001RW** Nitrous Oxide 4 Bar Input (94mm wide)
- T1002R** MA-7, 7 Bar Input (42mm wide)
- T1002RW** MA-7, 7 Bar Input (94mm wide)
- T1003R** MA-4, 4 Bar Input (42mm wide)
- T1003RW** MA-4, 4 Bar Input (94mm wide)
- T1004R** Vacuum (42mm wide)
- T1004RW** Vacuum (94mm wide)
- T1005R** Entonox® 4 Bar Input (42mm wide)
- T1005RW** Entonox® 4 Bar Input (94mm wide)

OVERVIEW

Oxylitre Holdings Limited offer the most comprehensive range of Pipeline and Hose Accessories. Our Pipeline and Hose Accessories cover all Medical Gases suitable. Our connections are suitable for use with all flexible medical hose. Our extensive range of gas specific pipeline connections and accessories cover both British and European Standard connections. A full and comprehensive look at our range can be found in our product index.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom
+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

