

PIPELINE ACCESSORIES

Oxygen 4 Bar Input (95mm wide)

PRODUCT CODES

T1000Y to T1005Y



SPECIFICATION

- Brass Chrome Plated Body
- NOT MRI Conditional
- British Standard Probe
- BS Indexed Coloured Self Sealing Valve

RANGE

- T1000Y** Oxygen 4 Bar Input (95mm wide)
- T1001Y** Nitrous Oxide 4 Bar Input (95mm wide)
- T1002Y** MA-7, 7 Bar Input (95mm wide)
- T1003Y** MA-4, 4 Bar Input (95mm wide)
- T1004Y** Vacuum (95mm wide)
- T1005Y** Entonox® 4 Bar Input (95mm wide)

OVERVIEW

Oxylitre Holdings Limited offer the most comprehensive range of Pipeline and Hose Accessories. Our Pipeline and Hose Accessories cover all Medical Gases suitable. Our connections are suitable for use with all flexible medical hose. Our extensive range of gas specific pipeline connections and accessories cover both British and European Standard connections. A full and comprehensive look at our range can be found in our product index.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom
+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

