


PIPELINE AND HOSE ACCESSORIES

N.I.S.T Couplings (Nut & Liner)

PRODUCT CODES

NI0020 to NI0025



SPECIFICATION

- Stainless Steel
- BS Indexed
- MRI Conditional
- 6.5mm O/D Tubing Tail

RANGE

NI0020 Oxygen

NI0021 Nitrous Oxide

NI0022 Entonox®

NI0023 Medical Air 4 Bar (MA-4)

NI0024 Medical Air 7 Bar (MA-7)

NI0025 Vacuum

OVERVIEW

Oxylitre Holdings Limited offer the most comprehensive range of Pipeline and Hose Accessories. Our Pipeline and Hose Accessories cover all Medical Gases suitable. Our connections are suitable for use with all flexible medical hose. Our extensive range of gas specific pipeline connections and accessories cover both British and European Standard connections. A full and comprehensive look at our range can be found in our product index.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

