


MEDICAL RAIL SYSTEMS AND FITTINGS

Universal Medical Rail System

PRODUCT CODES
8 codes



SPECIFICATION

	Extruded Aluminium
	Anodised Silver Finish
	MRI Conditional
	Fixes Directly to wall
	Complete with End Caps & Infill Strip
	Compatible with standard rail-mounted accessories
Lengths	500mm, 750mm, 1000mm, 1500mm, 2000mm, 2500mm, 3000mm, 6000mm
	Mounting holes are to be drilled upon installation

RANGE

R500	500mm
R750	750mm
R1000	1000mm
R1000-5	1500mm
R2000	2000mm
R2500	2500mm
R3000	3000mm
R6000	6000mm

ACCESSORIES

8 CODES

R1000A	End Cap Only (each)
R1000B	Rail Infill Strip Only (per meter)
R1302	Rail Clamp 44mm Only

OVERVIEW

Oxylitre Holdings Limited Rail Systems and Rail Fittings cover a comprehensive range of requirements. Our extruded heavy duty aluminium rail is suitable for use in various medical environments, and comes supplied with plastic end caps and infill. Anodised silver, our rail is available in numerous standard lengths, non-standard lengths are available on request. Our extensive range of Rail Fittings have been designed to create flexible solutions for the safe storage of Medical Sundries.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom
+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

