

BRITISH STANDARD REMOTE PROBES

British Standard Remote Probes (straight tube)

PRODUCT CODES

6 codes



SPECIFICATION

- Stainless Steel
- One Pieced Machined
- BS Gas Indexed
- MRI Conditional
- 6.5mm O/D straight Tubing Tail

RANGE

- 0032-S** Oxygen
- 0042-S** Nitrous Oxide
- 0052-S** Vacuum
- 0063-S** Medical Air 7 Bar (MA-7)
- 0064-S** Medical Air 4 Bar (MA-4)
- 00720-S** Entonox®

WEIGHTS AND DIMENSIONS

Dimensions

6.5mm OD Tube Tail to BS5682

OVERVIEW

Oxylitre Holdings Limited offer the most comprehensive range of Pipeline and Hose Accessories. Our Pipeline and Hose Accessories cover all Medical Gases suitable. Our connections are suitable for use with all flexible medical hose. Our extensive range of gas specific pipeline connections and accessories cover both British and European Standard connections. A full and comprehensive look at our range can be found in our product index.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom
+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • oxylitre.co.uk

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

