

BRITISH STANDARD SELF SEALING VALVES

British Standard Self Sealing Valves with 3/8 bsp Female Fitting

PRODUCT CODES

6 codes



SPECIFICATION

- Robust Moulded Construction
- Engraved Gas Identification
- BS Gas Indexed
- Colour Coded for Gas Service

RANGE

0085-38 Oxygen

0086-38 Nitrous Oxide

0087-38 Vacuum

0088-38 Medical Air 4 Bar (MA-4)

0089-38 Medical Air 7 Bar (MA-7)

0090-38 Entonox®

OVERVIEW

Oxylitre Holdings Limited offer the most comprehensive range of Pipeline and Hose Accessories. Our Pipeline and Hose Accessories cover all Medical Gases suitable. Our connections are suitable for use with all flexible medical hose. Our extensive range of gas specific pipeline connections and accessories cover both British and European Standard connections. A full and comprehensive look at our range can be found in our product index.

Morton House, Skerton Road, Old Trafford, Manchester M16 0WJ, United Kingdom

+44 (0)161 872 6322 • sales@oxylitre.co.uk • service@oxylitre.co.uk • [oxylitre.co.uk](https://www.oxylitre.co.uk)

Oxylitre Holdings Limited · Registered in the United Kingdom · Reg No. 11095954 · VAT No. GB 146 6892 29
MHRA Registration No. 7118 · EUDAMED SRN GB-MF-000021789 · © Oxylitre Holdings Limited 2026. All rights reserved.

